


# baebies® FINDER™



## Near-Patient Newborn and Pediatric Testing

Life-threatening conditions for newborns require multiple tests and large sample volumes for an accurate diagnosis. **FINDER** uses one drop of blood to rapidly perform critical tests.

### FINDER Features

- **3 analytes from a single test**
  - Glucose-6-Phosphate Dehydrogenase (G6PD)
  - Total Serum Bilirubin
  - Albumin
- **Small sample size**  
1 drop of whole blood (50µL)
- **Small footprint**  
12.4" L x 8.9" W x 10.6" H
- **Rapid time to result**  
<15 minutes
- **Simple workflow**  
No sample prep.



INSERT CARTRIDGE



ADD SAMPLE



RUN TEST



VIEW RESULTS

**For more information, contact us today!**

[info@baebies.com](mailto:info@baebies.com)

+1 (919) 891-0432

Durham, NC USA

[www.baebies.com/FINDER](http://www.baebies.com/FINDER)

everyone **deserves a healthy start**®

FINDER is not available at this time for sale or use in the USA or any other territory. The product has not been evaluated by the FDA or other regulatory authorities.  
Revision B

