

List of External Controls Qualified with FINDER® 1.5 Instrument and FINDER SARS-CoV-2 Test System

The FINDER 1.5 instrument and FINDER SARS-CoV-2 Test System, as part of lot release, use external controls manufactured by ZeptoMetrix LLC (878 Main Street, Buffalo, NY 14202 USA). These ready to use, inactivated full process controls are designed to evaluate performance of molecular tests. They are used by Baebies for validation of assays, training of laboratory personnel and to monitor assay-kit lot performance.

Within the Scope of this document:

Catalogue Name	Virus /Cell line	Target Concentration Copies/ mL
NATSARS(COV2)-ERC	SARS-COV-2 isolate USA-WA1/2020	50,000 Copies/ mL
NATSARS(COV2)-NEG	A-549	50,000 Cells/ mL

Qualified NATSARS(CoV2)-ERC Lots

Catalogue Name	Lots
NATSARS(CoV2)-ERC	325436
	325642
	325945
	327602

Qualified NATSARS(CoV2)-NEG Lots

Catalogue Name	Lots
NATSARS(CoV2)-NEG	324759
	324407
	325169
	326673

Please contact ZeptoMetrix LLC, 878 Main Street, Buffalo, NY 14202 USA • Tel (800) 274-5487 for additional information relating to the control materials.